



Lab Meeting





1. Why Create the Computational Model?




Traditional techniques

- Macroscopic release profiles
- Lack microscopic insights







The Carrier: Montmorillonite (MMT)



Act as a nano-sandwich


- a cost-effective, biocompatible excipient
- to improve the solubility of hydrophobic drugs.

Cyclosporine in Drug-Clay Complex





- An immunosuppressive agent
- Poor oral bioavailability(25%)
- Because of the absorption in the GI tract



Timeline of Drug Release



A. Deep Trapping



B. Edge Trapping



2. Preparation of Bioproduct


- Materials:**

Cyclosporine A (CsA): **0.25 g**

Montmorillonite (MMT): **1.25 g**

- Ratio:** 1:5 (w/w)

- Solvent:** 100 ml Ethanol (EtOH)



Using magnetic
stirrer for **24 hours**



- Equipment:** Rotary Evaporator
- Condition:** Water bath at **45°C**
- Process:** Evaporation under reduced pressure(**150 mbar**)


Step 1: Mixing

Step 2: Stirring

Step 3: Solvent Removal

XRD Analysis Results

X-Ray Diffraction (XRD) Characterization of Bioproducts




Cyclosporine + Clay


$d(001) = 12.04 \text{ \AA}$

Montmorillonite

$d(001) = 11.90 \text{ \AA}$




did not fully enter the interlayer space



Dissolution Test & HPLC Analysis

1. Dissolution Medium Preparation



Step A: Preparation of SIF (pH 6.8)



- Composition:**
 - KH₂PO₄: **3.402 g**
 - NaOH: **0.448 g**
 - Distilled Water: **0.5 L**
- Adjustment:** pH adjusted to **6.8**
(Simulated Intestinal Fluid).

Step B: Modification for Hydrophobic Drug


- Due to the low aqueous solubility of Cyclosporine, organic solvent was added to maintain **sink conditions**.
- Medium:** **60% SIF + 40% Ethanol (v/v)**.

Dissolution Test & HPLC Analysis

- Conditions:**
 - Temperature: **37 degrees**
 - Stirring Speed: **100 rpm**
- Reference (Vessels 1-3): 25 mg Pure Cyclosporine Capsules.** 🍯
- Test (Vessels 4-6): 0.15 g Bioproduct (Equivalent to 25 mg drug)*.** 🍯



**Volume: 1000 mL
Medium per vessel.**




- Sampling Points:** 10, 30, 60, 90, 120, 180, 240, 300, 360, 420, 480 min.
- Aliquot Volume:** **5 mL** withdrawn at each time point (replaced with fresh medium).
- Quantification:** Samples filtered and analyzed via **HPLC**.

Dissolution Test- Result

% liberado


—●— CYCLO only —●— CYCLO-MMT




Lower Release
-the lack of intercalation
(from XRD)

Large SD
-Heterogeneity in the
physical mixture
- More replicas are
required

The other Clay options



The structure of Sepiolite



The structure of Halloysite

Step 1: Mixing

Step 2: Stirring


Step 3: Solvent Removal

1. Medium
Preparation

2. Experimental
Setup

3. Sampling Protocol

The other Clay Sepiolite option- Dissolution Test



-immediate release, reaching
~**90%** within the first **10 minutes**.

-Reaches nearly **100%** equilibrium
quickly (within 30 mins)